



OCTANE

TITANIUM TECHNOLOGY:

A NEW LEVEL OF DURABILITY

SUB4



*Excluding wheels and cushion

OCTANE SUB4

RGK 30
1988 - 2018

THE WORLDS FIRST MADE TO MEASURE TITANIUM WHEELCHAIR GUARANTEED UNDER 4 KG*

The Octane Sub4 elevates the daily wheelchair to a whole new level of performance and engineering down to the finest detail.

We took a classic design and redefined it even further by combining Sub4 technology with a hand-crafted made to measure titanium grade 9 frame. Resulting, in an ultra-lightweight and rigid wheelchair that is perfectly designed to fit you and your lifestyle.

TITANIUM GRADE 9:

- STRENGTH TO WEIGHT RATIO
- SMOOTH RIDING
- NO CORROSION
- DURABILITY





Frame finish options: polished or brushed



Sub4 technology –
guaranteed under 4 kg*

KEY FEATURES

- Made to measure
- Less than 4 kg guaranteed (excluding wheels and cushion)
- Rigid and lightweight
- Titanium grade 9 frame
- Low profile castor arm design
- Airtech Sub4 technology
- Low profile cross braces
- Optimised single tube design
- Max user weight 110 kg



Low profile castor arm design



Single post axle

OCTANE

TECHNICAL DATA

SUB4

OCTANE SUB4 // Titanium Wheelchair

Max. user weight	110 kg	Total length	Custom-made
Frame version	Fully welded open frame	Total height	Custom-made
Frame inset	Custom-made	Max. pivot width	Custom-made
Frame angle	Custom-made	Max. safe slope	7° (for use of parking brakes)
Seat width	Custom-made	Front wheel	3" / 4" / 5"
Seat depth	Custom-made	Rear wheel	24" / 25"
Seat height front	Custom-made	Wheelchair weight	Less than 4 kg guaranteed
Seat height rear	Custom-made	Mass of heaviest part	2.1 kg
Backrest height	Custom-made	Removable parts	Rear wheels (for transport only)
Backrest angle	Custom-made	Colours	- Frame finish options: Brushed frame or polished frame - 5 accent colours (decal on frame), 6 backrest upholstery stitching colours, 9 Spinergy colour options
Centre of gravity (COG)	Custom-made	Intended use as seat in motor vehicle	Not released as seat in motor vehicles
Wheel camber	1° / 2° / 3°	Transport without user	Fix the wheelchair safe in the vehicle (e.g. car or plane)
Total width	SW + 240 mm	Intended use and environment	For people with limited mobility (children & adults); self-propelling or being pushed by an assistant, indoor and outdoor use.



Sunrise Medical Australia Pty Ltd.

6 Healey Circuit,
Huntingwood NSW 2148

Phone: (02) 9678 6600

Fax: (02) 9423 6969

Email: enquiries@sunrisemedical.com.au

www.SunriseMedical.com.au



For further information on the full specification, options and accessories, please refer to the order form. All information is subject to change without notification. Please consult Sunrise Medical with any queries you may have. If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com.au.

